




Elsam Dental Treatment Wet Lift Table Specifications (Back Views)



MOUNTING LOCATIONS

Hot and cold water supply can be located on the wall in the open area behind the table or on the wall next to the table using a shower control valve with a spray hose attached. Call to discuss options.



Drain can be located in this area of floor behind the table- (floor drain is preferred but wall drain is acceptable also.)

Dental Focus™ LLC
P.O. Box 454, Neshanic Station NJ 08853
e-mail: jmerritt@dentalfocus.biz www.dentalfocus.biz
Ph: (908) 230-1955