



ALL CLINICAL APPLICATIONS ARE EASY WITH THE ALL-IN-ONE AIR POLISHER

airN GO
easy

AIR-N-GO® SUPRA NOZZLE *(standard)*



- Suitable for any supra-gingival prophylactic treatment in combination with AIR-N-GO® "CLASSIC" Peppermint powder
- Faster than using prophylactic paste
- Treat fissures, grooves, and occlusal surfaces

AIR-N-GO® PERIO easy NOZZLE *(optional Perio easy kit)*



- "Shoehorn" shape providing easy sub-gingival insertion
- Disrupt biofilm below the gum line to treat periodontitis
- 5 seconds per face is all that is needed



AIR-N-GO® nozzles provide continuous, accurate and controlled spray. These nozzles offer power and efficiency for the removal of supra and subgingival biofilm.



AIR-N-GO® easy SUPRA 120°
 2 nozzles supplied as standard




AIR-N-GO® PERIO easy
 Nozzle supplied with the PERIO kit (optional)

Supragingival Powder

- Instant results
- Finishing after scaling
- Cleaning of interproximal areas, fissures and crevices
- Remove colored stains and bacterial plaque, even in areas that are difficult to access
- Gentle and efficient with limited bleeding during treatment

"CLASSIC"
 Sodium bicarbonate-based
 Average particle size 76µm
 FLAVORS: Peppermint




Particles with softened geometry (x1000)



Multiple Advantages with Unique Properties

- Plaque removal: during periodontal maintenance on exposed root surfaces without sacrificing substance
- Reduction of the bacterial load in the periodontal pockets
- Removal of the sub-gingival biofilm
- Treatment of periodontal disease

"PERIO"
 Glycine-based
 Average particle size 25µm
 FLAVOR : Neutral



AIR-N-GO® "PERIO" (x1000)
 Source: ICMCB CNRS (Institute of Condensed Matter Chemistry - French National Centre for Scientific Research) S.E.M. photos

