



P.O. Box 454, Neshanic Station, NJ 08853 | (908) 968-3553 | www.dentalfocus.biz

Sopix Minimum Configuration

- Operating System Windows 10
- Processor Intel Core i5
- RAM 6 GB
- Hard Drive 1 TB Hard Disk Drive
- Graphic Card Nvidia or ATI 2 GB Unshared Memory
- Screen Resolution 1280 x 1024 or More
- Ethernet Board 1Gbps
- Power Supply 200 Watts or Higher

Sopix Recommended Configuration

- Operating System Windows 10 Professional
- Processor Intel Core i7
- RAM 8 GB
- Hard Disk 250GB Solid State Drive
- Graphic Card Nvidia or ATI 2 GB Unshared Memory
- Screen Resolution 1920 x 1080 or More
- Ethernet Board 1 Gbps
- Power Supply..... 200 Watts or Higher

Important: Computer should be in view when shooting images. This allows you to see images as you work on cases and see mistakes instantly. We suggest wireless mouse and keyboard or a computer in reach when shooting. Walking across a room is not good work flow and can lead to damaged equipment.